

User Manual

Pneumonia Detection

Author: Vikas Ramaswamy
Matriculation No.: 7024164
Course of Studies: Masters in Business Intelligence & Data Analytics

Author: Adarsh Pal
Matriculation No.: 7024418
Course of Studies: Masters in Business Intelligence & Data Analytics

Author: Sangram Patil
Matriculation No.: 7023678
Course of Studies: Masters in Technical Management

First examiner: Prof. Dr. Elmar Wings
Submission date: July 3, 2023

Contents

1	Getting to know about Pneumonia Detection	1
2	Operation Of Pneumonia Detection	3
2.1	Software Setup	3
2.1.1	System Requirements	4
3	First Step in Pneumonia Detection	5
3.1	Interface for Pneumonia Detection:	5
3.2	Steps for Pneumonia Detection	5
4	Main Function	9
5	Help	11
5.1	Frequently Asked Questions	11
5.2	Troubleshooting	12
6	Functions of Pneumonia Detection Interface	15
7	Maintenance for Pneumonia Detection Web Application	17

1 Getting to know about Pneumonia Detection

What is Pneumonia Detection?

Pneumonia is a respiratory infection caused by bacteria or viruses; it affects many individuals, especially in developing and underdeveloped nations, where high levels of pollution, unhygienic living conditions, and overcrowding are relatively common, together with inadequate medical infrastructure. Pneumonia causes pleural effusion, a condition in which fluids fill the lung, causing respiratory difficulty. Early diagnosis of pneumonia is crucial to ensure curative treatment and increase survival rates. Chest X-ray imaging is the most frequently used method for diagnosing pneumonia. However, the examination of chest X-rays is a challenging task and is prone to subjective variability. In this study, we developed a computer-aided diagnosis system for automatic pneumonia detection using chest X-ray images.

Pneumonia detection involves identifying the presence of pneumonia, which is an infection or inflammation of the lungs. It typically includes a combination of clinical assessment, medical imaging (such as chest X-rays or CT scans), laboratory tests (such as blood tests and sputum cultures), and symptom evaluation. The goal is to identify characteristic symptoms, abnormal lung sounds, and imaging findings that indicate the presence of pneumonia. Additionally, the analysis of patient data, including demographics, medical history, and vital signs, helps healthcare professionals in making an accurate diagnosis and determining appropriate treatment options.

Pneumonia detection involves analyzing medical images, such as chest X-rays or CT scans, to identify signs of pneumonia. Here are some key features commonly used in Pneumonia Detection:

- **Image Preprocessing:** Prior to analysis, the medical images undergo preprocessing steps to enhance the relevant features and remove noise or artifacts. This may involve resizing, cropping, normalization, or applying filters to improve image quality.
- **Lung Segmentation:** In order to focus the analysis on the lung area, lung segmentation is performed to separate the lung region from the rest of the image. This step helps isolate potential pneumonia regions for further examination.
- **Chest X-ray:** Imaging plays a crucial role in pneumonia detection. Chest X-rays can help visualize the lungs and identify areas of infection or inflammation. Typical findings include infiltrates (consolidation or patchy opacities) in the affected lung regions.

1 Getting to know about Pneumonia Detection

- **Artificial Intelligence (AI) Algorithms:** The app can employ AI algorithms, such as machine learning or deep learning models, to aid in pneumonia detection. These algorithms can be trained on large datasets of pneumonia images to enhance accuracy and speed in identifying potential cases.
- **Accessibility and Security:** Accessible to ensure it can be used on various devices and browsers. Additionally, robust security measures should be implemented to protect patient data and comply with privacy regulations.

2 Operation Of Pneumonia Detection

2.1 Software Setup

To set up the software for pneumonia detection, you can follow these steps:

1. **Download and install Python:**

Ensure that Python is installed on your system. You can download the latest version of Python from the official Python website and follow the installation instructions specific to your operating system. Visit the official Python website (<https://www.python.org>) and download the latest version of Python for your operating system.

2. **Install PyCharm:**

Visit the JetBrains website (<https://www.jetbrains.com/pycharm>) and download PyCharm Community Edition, which is the free version. Install PyCharm by following the installation instructions specific to your operating system.

3. **Install Required Libraries:**

Use pip, the Python package manager, to install the necessary libraries and frameworks for pneumonia detection. Some essential libraries for image processing and machine learning in this context may include OpenCV, Pillow, TensorFlow, Keras.

4. **Set Up Flask Web Framework:**

Install the Flask web framework, which will be used to develop the web application. Use pip to install Flask

5. **Train Pneumonia Detection Model:**

Train your model using labeled data. If you train your model, ensure that the training data is properly prepared and augmented, and follow standard machine learning practices for training, validation, and evaluation.

6. **Build the Web Application:**

Develop the web application using Flask and the necessary HTML, CSS, and JavaScript files. Define routes, handle user requests, and integrate the pneumonia detection model into the application's backend.

2 Operation Of Pneumonia Detection

7. Test the Application:

Test the functionality of the web application to ensure that it correctly processes uploaded images, performs pneumonia detection, and provides accurate results. Use sample images or a test dataset to verify the performance of the application.

8. Deployment:

Choose a suitable platform or hosting service to deploy the web application. Popular options include cloud platforms like Google Cloud Platform (GCP)

2.1.1 System Requirements

- **Software requirements:**

Table 2.1: Software requirements

Software/Package	Open Source/License	Version
VS Code	Open Source	1.64.2
Anaconda-Navigator	Open Source	2.1.4
Conda	Open Source	4.11.0
Keras	Open Source	2.8.0
TensorFlow	Open Source	2.8.0

- **Hardware requirements:**

Table 2.2: Hardware requirements

Hardware Name	Description	Quantity
Computer Display	1080p (1920x1080) resolution or higher	1
Keyboard	Standard American/German language	1
Mouse	Standard 3 button wired or wireless	1
Laptop	Capable of running deep learning algorithms	1

3 First Step in Pneumonia Detection

3.1 Interface for Pneumonia Detection:




Figure 3.1: Interface Of Pneumonia Detection

3.2 Steps for Pneumonia Detection

Step 1: Access the Web Application:

Open a web browser and enter the URL or click on the provided link to access the pneumonia detection web application.

Step 2: Upload Chest X-ray Image:

Click on the "Choose File" button to open a file selection dialog box.

Step 3: Select the Chest X-ray Image:

In the file selection dialog box, navigate to the location where the chest X-ray image is stored on your device. Select the desired image file and click on the "Open".

3 First Step in Pneumonia Detection




Figure 3.2: User Interface




Figure 3.3: Choose File Option




Figure 3.4: Selection Of File

Step 4: Image Processing and Analysis:

After selecting the image click on the “Submit”. After the image is uploaded, the web application will start processing. This process may take a few seconds to complete.




Figure 3.5: Submit Option

Step 5: View Pneumonia Prediction:

Once the analysis is finished, the web application will present the pneumonia prediction result to you as shown in Figure.




Figure 3.6: Result Display

Step 6: Interpret the Result:

There are two possibility that is “Person in Healthy” or “Person has Pneumonia”

- **Positive Result:**

For this, message is displayed as “Person has Pneumonia”.As shown in fig 3.7

- **Negative Result:**

For this the, message is displayed as “Person in Healthy”. As shown in fig. 3.8

3 First Step in Pneumonia Detection




Figure 3.7: Positive Result




Figure 3.8: Negative Result

4 Main Function

Main Function

Pneumonia causes pleural effusion, a condition in which fluids fill the lung, causing respiratory difficulty. Early detection of pneumonia is essential for ensuring curative care and boosting survival rates. The approach most usually used to diagnose pneumonia is chest X-ray imaging. However, examining chest X-rays is a difficult task that is vulnerable to subjectivity. The major purpose of the Web app for pneumonia detection is to give consumers a place to upload chest X-ray images and get assessments of whether pneumonia is present or not. The Pneumonia detection web application analyses the provided images and produces prediction results using an image processing pipeline and a trained deep learning model.

An outline of the primary actions that make up the Pneumonia Detection operation is provided below




Figure 4.1: User Road Map

1. Image Upload:

Through a simple user interface, the web application enables users to upload their chest X-ray images. To choose the image file from their local device, users often click on a particular area or button.

2. Image Preprocessing:

Segmentation separates the lung region for study while preprocessing techniques improve the image quality. A classification algorithm, such as machine learning or deep learning models, determines the presence or absence of pneumonia by extracting the necessary information from the segmented images. Additionally, it sends notification to fix the uploading image, indicating the precise issue, if the image is not in the suitable format, such as size, blurry image, or incorrect file format.

4 Main Function

3. Prediction Generation:

After being preprocessed, the image is then entered into a deep learning model that has been trained particularly to detect pneumonia. The model evaluates the image and makes a determination regarding the presence or absence of pneumonia. Predictions can be made in the form of messages such as "You have pneumonia or you don't."

4. Result Display:

The user interface of the web app shows the user the prediction result. The result is typically presented in the form for message. Stating that the use have pneumonia or not




Figure 4.2: Result

5 Help

5.1 Frequently Asked Questions

1. How do I upload an image?

To upload an image, click on the "Upload" button or drag and drop the image file into the designated area. Make sure the image file is in a compatible format, such as JPEG or PNG.

2. What should I do if my image is not uploading?

Ensure that the image file size does not exceed the specified limit. If the problem persists, try using a different web browser or check your internet connection.

3. Can I upload multiple images at once?

Currently, the web app supports the upload of one image at a time. However, you can repeat the process for multiple images individually.

4. How long does it take to get the prediction results?

The prediction results are generated within a few seconds after the image is uploaded. The processing time may vary depending on the size of the image and the server's processing capacity.

5. How accurate are the prediction results?

Our model has been trained on a large dataset and has shown high accuracy in pneumonia detection. However, it's important to note that no diagnostic tool is 100% accurate, and consulting a medical professional is always recommended for accurate diagnosis.

6. How do I understand the outcome of the prediction?

If the prediction result indicates "Pneumonia Detected," it means that the model has detected pneumonia in the uploaded image. If the result indicates "No Pneumonia Detected," it means that pneumonia is not detected. However, it's important to consult with a medical professional for an accurate diagnosis and further medical advice.

5.2 Troubleshooting

Troubleshooting is an essential aspect of maintaining an image processing web application for pneumonia detection. Here are some common issues that users may encounter and possible solutions to address them:

1. Image Upload Failure:

- Verify that the uploaded image file is in a supported format (e.g., JPEG, PNG) and meets the size restrictions set by the web app.
- Check the internet connection to ensure a stable and reliable connection.
- Clear the browser cache and try uploading the image again.
- If the issue persists, provide an error message to the user with instructions on how to resolve the problem or contact support for further assistance.

2. Slow Processing Time:

- Check the server resources and ensure they are sufficient to handle the image processing workload.
- Optimize the image processing pipeline by implementing efficient algorithms or techniques to speed up the analysis.
- Inform the user about the estimated processing time to set expectations and reduce frustration.

3. Inaccurate or Unreliable Results:

- Ensure that the deep learning model used for pneumonia detection is trained on a diverse, updated and high-quality dataset.
- Regularly update the model by incorporating new data and retraining it to improve accuracy.
- Implement cross-validation techniques to assess the model's performance and identify areas for improvement.
- Consider including a confidence score or probability value with the prediction results to indicate the level of certainty.

4. Compatibility Issues:

- Verify that the web application is compatible with different browsers and operating systems. Test the application on various platforms to identify and resolve any compatibility issues.
- Check for any dependency conflicts and ensure that the required libraries and frameworks are up to date.
- Provide clear instructions to users regarding the recommended browsers or system requirements for optimal performance.

5. Error Handling and Reporting:

- Implement comprehensive error handling mechanisms to catch and report any unexpected errors or exceptions.
- Log detailed error messages to facilitate debugging and troubleshooting.
- Display user-friendly error messages to guide users and provide instructions on how to resolve common issues.
- Set up a system to collect user feedback and error reports, allowing you to identify recurring issues and address them promptly.

6. User Support :

- Provide a dedicated support channel, such as an email address or online chat, for users to reach out with their concerns or issues.
- Maintain up-to-date documentation and user guides that cover common troubleshooting scenarios and provide step-by-step instructions for resolving them.
- Communicate known issues or limitations clearly to users, along with any workarounds or alternative solutions.

6 Functions of Pneumonia Detection Interface

Functions of Pneumonia Detection Interface

1. Choose File:

To start the analysis for the detection of pneumonia, click "Choose File" to upload your chest X-ray image. It also displays information about the file's selection state as shown in the figure.




Figure 6.1: Coose File Function

2. Submit:

The processing and analysis of the submitted image for the identification of pneumonia are started using the "Submit" feature. The user can view the results of the prediction by clicking the "Submit" button, which instructs the web app to process the image according to the trained model.




Figure 6.2: Submit Function

3. Go Back Functions:

By clicking the "Go back" button, the user can navigate back to a previous screen. Where user can upload the other photo for analysis.




Figure 6.3: Go back Function

4. Result Visualization:

The result visualization function is responsible for displaying the prediction results to the user in a clear and informative manner. After the image processing and analysis are completed, the result visualization function presents the prediction outcome, indicating whether pneumonia is detected or not.

- Image : In this the uploaded image is showed
- Message: Message indicated the prediction result like “Person is Healthy” or “Person has pneumonia”




Figure 6.4: Result Visualization

By providing accurate and relevant information for each input parameter, you enable the model to generate more precise sales predictions tailored to your specific requirements and context.

7 Maintenance for Pneumonia Detection Web Application

As the online realm is ever-evolving, the Pneumonia Detection Project demands periodic maintenance and frequent upgrades to function optimally. An accurate and efficient operation of a web app for pneumonia Detection depends on routine upkeep and upgrades. Here are some important factors to be consider for maintaining a Pneumonia Detection web app.

1. **Model Updates:**

Maintain a close eye on developments in algorithms and research for the Pneumonia Detection in order to implement improvements and raise prediction accuracy. To keep the model performing well, periodically retrain it using new data.

2. **Dataset Updates:**

Ensure the dataset used for training the model is regularly updated. Incorporate new data that represents a diverse range of pneumonia cases to improve the model's ability to detect different variations of pneumonia accurately.

3. **Bug Fixes and Performance Optimization:**

Continuously monitor and address any bugs or issues reported by users. Conduct regular testing and performance optimization to enhance the web app's speed, efficiency, and user experience.

4. **Security Updates:**

Stay up to date with the latest security protocols and best practices to safeguard user data and protect against potential security vulnerabilities. Maintain the web app's safety measure by implementing the essential security protections and upgrades.

5. **Compatibility and Dependency Updates:**

Keep track of any changes or updates to the dependencies, libraries, or frameworks used in the web app. Regularly update these components to ensure compatibility with the latest versions and maintain a stable environment.

6. **Backup and Recovery:**

Implement a reliable backup plan to protect user data and ensure its recoverability in case of any unforeseen events or data loss. Regularly back up the relevant data and establish recovery procedures.

7 Maintenance for Pneumonia Detection Web Application

7. User Feedback and Support:

Maintain a support system to address user queries and provide assistance as needed related to Pneumonia Detection. Incorporate user feedback to enhance the web app's usability and features.

8. Documentation and User Guides:

Keep the documentation and user guides up to date to reflect any changes or updates made to Pneumonia Detection Web App . Ensure that users have access to accurate and relevant information on how to use the app effectively.